

DEFIBTECH ARM XR 7-YEAR WARRANTY

The newest addition to the Defibtech product family is the **ARM XR Automated Chest Compression (ACC) device**, designed for professional rescuers. ARM XR ACC helps emergency personnel deliver continuous, high-quality CPR for patients in cardiac arrest. Easy to deploy and use, ARM XR allows first responders to transition from manual to mechanical CPR in seconds and maintain continuous compressions over long durations. Count on ARM XR to deliver consistent, high-quality CPR from the rig to the hospital.



PURCHASE COMES PATIENT-READY WITH:

- Carry Case (D-shaped)
- Stabilization Strap
- Wrist Straps
- USB Cable
- Three (3) Suction Cups
- Backboard
- Frame Assembly
- Battery Pack
- Power Module/Charge

SPECS/PRODUCT DIFFERENTIATORS

- **Active Chest Recoil** - Providing up to 0.6 in (1.5 cm) of recoil for effective chest wall expansion
- **Consistent Resuscitation Compressions** - Delivering reliable, uninterrupted compressions
- **Integrated Frame and Backboard** - Enabling faster patient application with snap-and-release functionality
- **Long Battery Life** - Operating on battery or AC power
- **Compatibility with Many Patient Sizes** - Accommodating chest widths up to 17.5 in (44.4 cm) and chest heights from 7.4 to 12.7 in (18.8 to 32.3 cm)
- **Lightweight and Portable Design** - Allowing for easy handling with a max weight of 16.4 lbs (7.5 kg) including the battery pack.

WITH ARM XR, EMS TEAMS CAN:

- Run a cardiac arrest code with smaller teams and reduced fatigue/stress
- Treat large and small patients with automated compression depth adjustment
- Extend budgets with this cost-effective, rugged mechanical CPR system



DEFIBTECH ARM XR 7-YEAR WARRANTY

Every ARM XR comes with **MME's exclusive 7-year warranty**. Warranty covers parts and labor for all repairs needed due to manufacturer defect or faulty workmanship.

Warranty excludes repairs needed due to device abuse, misuse, or acts of God. Device must have Preventative Maintenance (PM) performed by ReNew Biomedical and ReNew Biomedical only. Service, Repairs and PMs performed at ReNew Biomedical Depot only. Warranty is void if equipment is diagnosed, altered, serviced, or repaired by any party other than ReNew Biomedical. Warranty does not apply to back plate, side supports, batteries, charging cables, cases, or other accessories. Warranty becomes null and void if annual PM is not performed within 60 days of PM due date.