

DEFIBTECH ARM XR

The newest addition to the Defibtech product family is the **ARM XR Automated Chest Compression (ACC) device**, designed for professional rescuers. ARM XR ACC helps emergency personnel deliver continuous, high-quality CPR for patients in cardiac arrest. Easy to deploy and use, ARM XR allows first responders to transition from manual to mechanical CPR in seconds and maintain continuous compressions over long durations. Count on ARM XR to deliver consistent, high-quality CPR from the rig to the hospital.



RENTAL COMES PATIENT-READY WITH:

- Carry Case (D-shaped)
- Stabilization Strap
- Wrist Straps
- USB Cable
- Three (3) Suction Cups
- Backboard
- Frame Assembly
- Battery Pack
- Power Module/Charge

SPECS/PRODUCT DIFFERENTIATORS

- **Active Chest Recoil** - Providing up to 0.6 in (1.5 cm) of recoil for effective chest wall expansion
- **Consistent Resuscitation Compressions** - Delivering reliable, uninterrupted compressions
- **Integrated Frame and Backboard** - Enabling faster patient application with snap-and-release functionality
- **Long Battery Life** - Operating on battery or AC power
- **Compatibility with Many Patient Sizes** - Accommodating chest widths up to 17.5 in (44.4 cm) and chest heights from 7.4 to 12.7 in (18.8 to 32.3 cm)
- **Lightweight and Portable Design** - Allowing for easy handling with a max weight of 16.4 lbs (7.5 kg) including the battery pack.

WITH ARM XR, EMS TEAMS CAN:

- Run a cardiac arrest code with smaller teams and reduced fatigue/stress
- Treat large and small patients with automated compression depth adjustment
- Extend budgets with this cost-effective, rugged mechanical CPR system

ADDITIONAL RENTALS INFORMATION:

No minimum or maximum rental length

Annual preventative maintenance included with 12 month contract

Rental units come patient-ready